



LED.B

Curing Light

Classic and Enduring

Wireless.



Time setting: 5s, 10s, 15s, 20s.



Large capacity battery.
A full charge can be used for more than 500 times continuously under 10s working time mode.

Low standby power consumption with seventy days' standby time.



Constant light intensity.
The solidification effect is not affected by the consumption of remaining power.



technical parameter

Power input: AC100V-240V 50Hz/60Hz
Light output: 1000mW/cm² -1200mW/cm²
Net weight: 145g
Dimensions: 260mmx31mmx34mm